



## Professional Office Neckloop System™



### The ClearSounds Professional Office Neckloop System™

*Bringing new levels of communications capabilities to the office for people with a hearing loss.*

The system includes the powerful and advanced designs included in:

- The popular Plantronics M22 Headset Base Amplifier providing full range frequency control, and a simple but elegant Compatibility Selector that virtually guarantees compatibility with any business telephone
- The ClearSounds CLA7v2 Amplified Handsfree Neckloop providing clear binaural or monaural amplification for t-coil equipped hearing aids and cochlear implants and connectivity cable options to most handsfree or audio devices using 2.5mm or 3.5mm jacks.
- The ClearSounds M22 Connectivity Cable that provides modular connectivity between the M22 Base and the CLA7v2 Neckloop.
- The handset adapter that allows use of either the existing handset or the CLA7v2 Neckloop with the M22 Base.

For people with hearing loss and who have t-coil equipped hearing aids or cochlear processors, the Professional Office Neckloop System™ provides unparalleled sound management capabilities. For binaural device users, the ability to hear thru both ears can improve speech discrimination by up to 80%.

The noise cancelling microphone used in the CLA7v2 minimizes ambient noise and maximizes the ability of the person on the other side of the call to hear you as clearly as possible.

The System is great for improving employment opportunities for people with hearing loss. With the wide-spread use of digital telephone systems and now VOIP systems, office telephone communications has been a significant problem for people with a hearing loss. The ClearSounds Professional Amplified Office Neckloop System™ provides the ability for accommodating most hearing loss requirements and enhancing communication capabilities.

Distributed by:



#### Sonic Innovations

1/10 Dorsey Street, Milton, QLD 4064

PO Box 1869, Toowong DC, QLD 4066

Free call: 1800 639 263 F 07 3368 1044

E enquiry@sonici.com.au W www.sonici.com.au