

BILL TO ACCOUNT # _____	SHIP TO ACCOUNT # _____	PATIENT INFORMATION:
NAME: _____	NAME: _____	FIRST NAME: _____
ADDRESS: _____	ADDRESS: _____	LAST NAME: _____
PHONE: _____	PHONE: _____	CLINICIAN CONTACT: _____
	SAME AS PER BILL TO ADDRESS: <input type="checkbox"/>	PURCHASE ORDER: _____
		DATE REQUIRED: _____

NEW ORDER: REPAIR: REMAKE: ORIGINAL SERIAL NUMBER:

AUDIOGRAM		Hz	250	500	750	1K	1.5K	2K	3K	4K	6K	8K
RIGHT	AC											
	BC											
LEFT	AC											
	BC											

Insert Phones Headphones

BTE /miniRITE Colour

Celebrate 100 miniBTE
 Celebrate 80 miniRITE
 BTEP

Colour	R	L
Beige		
Taupe		
Brown		
Light Grey		
DarkGrey		
Black		

CUSTOM INSTRUMENTS	RIGHT						LEFT							
	IIC*	CIC	CICP	ITCD	ITC PDW	ITED	ITE PDW	IIC*	CIC	CICP	ITCD	ITC PDW	ITED	ITE PDW
<input type="radio"/> Celebrate ¹⁰⁰														
<input type="radio"/> Celebrate ⁸⁰														
Directionality				S	S	S	S				S	S	S	S
Program Button														
Volume Control														
Telecoil				S	S	S	S				S	S	S	S
Auto-tel				S	S	S	S				S	S	S	S
Wireless						S	S						S	S

NOTE: Please tick the options you require, if left blank the standard options will apply
Key: * Celebrate100 IIC only

CELEBRATE THIN TUBES, RECEIVERS, DOMES RIGHT LEFT BOTH

Dome Options	miniBTE	6mm	8mm	10mm	12mm	Tube Length	Tube Diameter
Open Dome						<input type="radio"/> Zero	<input type="radio"/> 0.9 (Std)
Bass, Double Vent						<input type="radio"/> One	<input type="radio"/> 1.3
Bass, Single Vent						<input type="radio"/> Two	
Power Dome						<input type="radio"/> Three	

Receiver Length	<input type="radio"/> One	<input type="radio"/> Two	<input type="radio"/> Three	<input type="radio"/> Four
Speaker Type	<input type="radio"/> 60	<input type="radio"/> 85	<input type="radio"/> 100	

WIRELESS OPTIONS:

Sound Gate	Sound Gate Mic *	* requires SoundGate 3 and above
TV Adaptor	Phone Adaptor	
RC-N Remote		

VENT & CANAL OPTIONS

Vent Type	Vent Size	Canal Length
Trench (Size = Open)	Extra Large (3.0mm)	As Marked on Impression
Collection	Large (2.4mm)	Deep - at 2nd bend
IROS	Medium (1.4mm)	Long-before 2nd bend
Straight	Small (1.0mm)	Medium-at 1st bend
Pressure (Size = Small)		Short-before 1st bend
No Vent		

FACEPLATE/SHELLCOLOR:

Beige/Beige Mid Brown/Clear Black (IIC only)
 Taupe/Clear Dark Brown/Clear

OTHER OPTIONS:

Canal Lock Soft Coat Name on Shell Removal Cord (Std. on CIC, CICP) Removal Notch

SPECIAL INSTRUCTIONS:

EARMOULD OPTIONS

Earmould Style	Material			Vent			Tube		
	R	L		R	L		R	L	
Micromould V1*						1 mm			Thin Tube 0.9mm
Sock (sleeve) mould*						1.4 mm			Thin Tube 1.3mm
Canal						2 mm			Standard BTE Tube
Canal Lock						2.4 mm			* NOT available with standard BTE tube option
Half Shell						3.0 mm			
Full Shell (carved)									
Skeleton									
LITE tip									
P100 Power Mould									

Canal Length	Options	Earmould Colour
<input checked="" type="radio"/> Factory Select	<input type="radio"/> As Marked	<input checked="" type="radio"/> Clear
<input type="radio"/> Deep	<input type="radio"/> Medium	<input type="radio"/> Beige
<input type="radio"/> Long	<input type="radio"/> Short	

SHIPPING:

Standard VIP - 2 Days - \$60 / Instrument