

The Pocket Listener

Hands-free personal assistive listening device for use straight out of the box

The Pocket Listener is a powerful, lightweight Assistive Listening Device used together with Telecoil-equipped hearing instruments. It provides users with greatly enhanced sound level and clarity from nearby sound sources or when connected directly to audio-visual equipment.

Designed for easy use straight out of the box with either single or binaural hearing instruments, The Pocket Listener amplifies the available sound level and noticeably improves the sound's clarity for greater hearing pleasure, even in noisy environments.

Easy to use – First, users select the Telecoil setting on their hearing instrument(s), then simply place The Pocket Listener around their neck and switch it on. The volume can be adjusted to a comfortable level by using the adjustment wheel on the side of the unit. Users will immediately notice an improvement in sound quality, whether conversing with friends and family, out shopping, or listening to a performance or speaker up to 10 meters away.

For enhanced listening pleasure to audio devices, users can also connect The Pocket Listener directly to their radio, stereo/HiFi, computer, or MP3 player with the jack-to-jack cable conveniently provided with the

The Pocket Listener – at a glance:

- Compact, lightweight hands-free design
- Worn around the neck, above or beneath clothing
- Easy use with single or binaural hearing instruments (with a Telecoil setting)
- Enhances sound level and clarity within a range of 5-10 meters
- Plugs directly into audio equipment (jack-to-jack cable provided)
- Rechargeable battery provides 8 hours of continuous use
- One year guarantee



unit. The Pocket Listener can also be connected to a television, DVD and video players, or Cable TV boxes using the optional TV Scart Adaptor Cable, which is sold separately.



'LVWULEXWHG E\

SONIC
Everyday Sounds Better

1/10 Dorsey Street, Milton, QLD 4064
PO Box 1869, Toowong DC, QLD 4066
Free call: 1800 639 263 F 07 3368 1044
E enquiry@sonici.com.au
W www.sonici.com.au

POCKET LISTENOR TECHNICAL INFORMATION

The Pocket Listenor is a powerful lightweight Assistive Listening Device used together with a Tele-coil equipped hearing instruments. It provides users with greatly enhanced sound level and clarity from nearby sound sources or when connected directly to audio visual equipment.

Physical Data

L x W x H: 65 x 34 x 25mm
Weight: 55g

Environmental

Case Material	Black moulded ABS
Water/dust intrusion rating	IP54
Drop test	100mm any axis
Operating temperature	0C – 70C
Storage temperature	-20C - +85C
Humidity	95% - non condensing

Electronic

Acoustic Sensitivity	65dB to 98dB SPL @ 1KHz
Distortion (THD) @ 1KHz	5.4%
S/N ratio	40dB
Neck loop output current	1A RMS
Supply Current Max vol	150mA
Supply Current Min vol	10mA
Supply Current Standby	50uA

Audio bandwidth	100Hz to 8KHz
Batteries	Ni MH 2-cell battery pack
Charger unit Input voltage:	100-240V 50Hz ac
Charge cycle time	8 Hours
80% charge time	2 Hour

Standards Compliance

C Council Directive 2004/108/EC(EMC DIRECTIVE)
FCC complies with Part 15 of the FCC Rules Operation

Line input Sensitivity	4mV RMS @ 50 kOhms
Input Mix	External input mutes mic
Battery Indicator Light	Light fades when voltage drops, indicating need to recharge.

