

Guymark VRA-002 Toy Box VRA System



A two cabinet toy box system for use in paediatric audiology with moving toys, set behind grey windows that “appear” when illuminated and can be independently animated.

Standard Features

- Illumination with or without animation
- Flashing or steady illumination without animation
- Mains operation with low voltage toys & control
- Supplied complete with cables and handset

Technical Specification

Cabinet size - 30cm x 30cm x 64cm
Power supply - 240v 50Hz or 110v 50-60Hz
Animation - 6v DC via 4 pin DIN

Toy view window size - 27cm x 43cm
Illumination - 6v DC via 4 pin DIN

Distributed by

One-stop-shop for all your audiological needs
www.sonici.com.au Free call 1800 639 263
Sonic Innovations Pty Ltd Unit 1, 10 Dorsey Street Milton Qld 4064
P (07) 3368 4044 F (07) 3368 1044 E enquiry@sonici.com.au

 **SONIC**
Everyday Sounds Better