

AVANT™ A2D+ Audiometer

Compact PC-Based Diagnostic Audiometer

The AVANT A2D+ Is A Powerful Dual Channel Audiometer That Features Integrated Quick SIN™ Testing Capability




Actual Size:
5" x 5" x 1.25"

MedRx®

www.medrx-usa.com or Call Toll Free: 888-392-1234

AVANT A²D⁺ Audiometer

Dual Channel Air, Bone & Speech Audiometer



The AVANT™ A²D⁺ Diagnostic Audiometer is a compact dual channel air, bone and word recognition system

- **Integrated Quick SIN™ Testing Capability**
- **Dual Channel**
- **Air, Bone & Speech Functions**
- **PC Based Via USB Connection**
- **Insert Earphones or TDH 39 Headset**
- **Patient Response Switch**
- **Ergonomic Design Available In Black Or White**
- **NOAH™, Sycle.Net & TIMS® Compatible**

MedRx

www.medrx-usa.com

1200 Starkey Rd. #105,

Largo, FL 33771

Toll Free 888-392-1234

727-584-9600 • Fax: 727-584-9602

A2D-S-ADH0107-1

Distributed by:

 **SONIC**
Everyday Sounds Better

1/10 Dorsey Street, Milton, QLD 4064
PO Box 1869, Toowong DC, QLD 4066
Free call: 1800 639 263 **F** 07 3368 1044
E enquiry@sonici.com.au **W** www.sonici.com.au