

Providing Audiological Solutions

Video Otoscope

Otoscope Probe, High Resolution Camera & Battery Operated LED Light Source

The light output of the LED produces virtually no heat or light outside the visible range. The LED light output is in a color range which can be easily white balanced by the MedRx video camera to render true color video otoscopy.

The remarkably clear resolution is particularly helpful in contrasting a clean ear canal and normal eardrum with an ear canal and eardrum that has excessive cerumen, foreign bodies, external otitis, PE tube status, cholesteatoma, TM perforations or other abnormalities.

The LED light source interfaces with the MedRx 3mm probe and miniature video camera which requires no focusing.



OS-S-ADH038-4



MedRx INC™

www.medrx-usa.com

1200 Starkey Rd. #105,

Largo, FL 33771

Toll Free 888-392-1234

727-584-9600 • Fax: 727-584-9602

MedRx's Authorized Representative in Europe
(Regulatory affairs only) Emergo Europe • Molenstraat 15
2513 BH The Hague • The Netherlands
Tel: (31) (0) 70 345-8570 • Fax: (31) (0) 70 346-7299

Video Companion

Bringing Video Otoscopy To Your Computer



The MedRx Video Companion is designed to easily interface with a laptop or desktop computer via a USB 2.0 video capture cable. This enables the clinician the ability to capture, view and store digital images in our NOAH™ compatible Video Otoscopy Software.

The primary components of the video otoscope are the deluxe probe, miniature video camera and the built in light source. The deluxe probe has an aluminum body and stainless steel conical tip that tapers down to 3mm. The miniature video camera is light-weight, streamlined, has a precision "white balance" built in and requires no focusing.

System Includes:

- Deluxe 3mm Otoscope Probe
- Compact Video Camera
- 1.25 Watt LED Light Source
- Battery Powered Light Source Includes Spare Batteries & Charger
- Focus Always Feature
- USB 2.0 Video Capture Cable
- Video Otoscope Software
- NOAH™ & TIMs® Compatible
- All Necessary Cables & Attachments
- 1 Year Warranty

Distributed by:



1/10 Dorsey Street, Milton, QLD 4064
PO Box 1869, Toowong DC, QLD 4066
Free call: 1800 639 263 F 07 3368 1044
E enquiry@sonici.com.au
W www.sonici.com.au