



## Are Your Directional Fittings on Thin Ice?

Directional microphones can provide hearing instrument users with great benefits if they are functioning properly. Now, in addition to comprehensive ANSI/IEC standards testing, the new AVANT Polar HIT™ allows hearing care professionals to verify the proper functionality of directional microphones prior to the initial patient fitting and after any future servicing. This advanced PC-based test chamber has a remarkably small footprint (11" w x 8.6" d x 4.8" h), yet is powerful enough to meet the needs of the most sophisticated practitioner.

- Unique rotating test pad for generating directional Polar Plots
- ANSI and IEC Hearing Instrument Testing
- Simple, one-cable USB connectivity
- Superior acoustic isolation for fast accurate testing



# MedRx INC™

Visit [www.medrx-usa.com](http://www.medrx-usa.com) or call toll free 888-392-1234 to learn more

If present climatic trends continue, Polar Bears could be extinct by the end of this century. With every purchase of an AVANT Polar HIT, MedRx will make a donation in your name for Polar Bear preservation.

Distributed by:

 **SONIC**  
Everyday Sounds Better

1/10 Dorsey Street, Milton, QLD 4064  
PO Box 1869, Toowong DC, QLD 4066  
Free call: 1800 639 263 F 07 3368 1044  
E [enquiry@sonici.com.au](mailto:enquiry@sonici.com.au)  
W [www.sonici.com.au](http://www.sonici.com.au)