

Clinical Audiometry **REDEFINED**



**Clinical Audiometry Has Reached
A New Level Of Performance
With The PC-Based AVANT
Stealth Audiometer Utilizing
Sound Card Technology**

MedRx[®]

www.medrx-usa.com or Call Toll Free: 888-392-1234

AVANT Stealth Audiometer

Clinical Audiometry Has Reached A New Level Of Performance



The AVANT™ Stealth Audiometer utilizes a single plug-and-play USB connection through the booth for trouble free connection to the computer

Cert #507707 CE
AD-S-ADH093-1

- Dual Channel
- PC Interface
- Internal Amplifiers
- Separate Transducer Outputs
- Compact Design:
12" x 11" x 2.25"
- NOAH™ & TIMS® Compatible
- High Frequency Option Available

Distributed by:



1/10 Dorsey Street, Milton, QLD 4064
PO Box 1869, Toowong DC, QLD 4066
Free call: 1800 639 263 F 07 3368 1044
E enquiry@sonici.com.au
W www.sonici.com.au