

OTOWAVE EARTIPS INFORMATION SHEET

The Otowave 102 is now supplied with the Type 1 eartips as shown below in place of the earlier Type 2 eartips.



Type 1



Type 2

Both sets of eartips may be used with the Otowave 102 but they ***must*** be used with the appropriate probe tip.

- **Type 1 uses the clear colourless probe tip**
- **Type 2 uses the opaque white probe tip**

| TYPE 1 | | | TYPE 2 | | |
|------------|---|--------|------------|---|--------|
| Stock Code | Size | Colour | Stock Code | Size | Colour |
| E9747 | 7mm | Blue | E9462 | 6mm | Grey |
| E9748 | 8mm | Yellow | E9463 | 7mm | Orange |
| E9749 | 9mm | Green | E9464 | 8mm | Green |
| E9750 | 10mm | Red | E9465 | 9.5mm | Blue |
| E9751 | 11mm | Blue | E9467 | 11mm | Grey |
| E9752 | 12mm | Yellow | E9468 | 12.5mm | Orange |
| E9753 | 13mm | Green | E9469 | 14mm | Green |
| E9754 | 14mm | Red | E9470 | 16mm | Blue |
| E9755 | 15mm | Blue | E9471 | 18mm | Grey |
| E9756 | 19mm | Yellow | | | |
| PROBE TIP | | | PROBE TIP | | |
| E9760 |  | Clear | E9761 |  | White |

N.B The colours referred to above are for illustrative purposes only and maybe subject to change.