

Quattro Bluetooth Amplified Accessory

Clear Features:

- 4 unique applications: Mobile, Home, Office, and Entertainment
- Accommodates T-coil and non-T-coil users
- Loop receiver requires T-coil(s)
- Headset jack accommodates audio headset for non-T-coil use (For best results, use the ClearSounds SmartSound earbuds)
- 30 dB adjustable amplification
- Compatible with standard Bluetooth Version 2.0 + EDR devices
- Profiles supported: Hands-Free Profile (HFP), Headset Profile (HSP) and Advanced Audio Distribution Profile (A2DP)

Key Features:

- Maximum range: 30 feet
- Battery performance: up to 8 hours talk time and up to 5 days standby
- Reduces or eliminates interference between cellular phones and hearing aid(s) or cochlear processor(s)
- Easy one-touch answering for phone calls
- Compatible with voice dialing (if supported by your cell phone provider)
- Multi-function key for headset volume control and "next song" selection during audio streaming

Additional Features

- Excellent reception and transmission of audio signals
- User-friendly controls & cables
- No interference from background noise
- T-coil compatible
- Binaural design
- Push-button slider to cinch neckloop for closer microphone placement
- Trouble-free operation
- Offers a 30 dB max gain
- Volume control for easy adjustment and multiple levels of amplification
- Built-in microphone
- Quick release feature
- Battery powered with a lithium ion rechargeable battery
- Charging options: AC adapter, USB cable and car charger (included)
- Size: 1.25" x 2.5" x .75"
- Color: Black
- 1-year warranty
- Item #: CS-QUATTRO



Compatibility

- Any standard Bluetooth enabled mobile phone
- Any standard Bluetooth enabled computer
- Any standard Bluetooth enabled cordless phone
- Any radio device that has a standard Bluetooth transmitter
- Certain iPods or MP3 players (using a standard Bluetooth transmitter)
- Other audio devices with the pairing device using v1.2 or v2.1

Quattro with ClearSounds RS062 SmartSound earbuds



Charging Accessories



Controls Background Noise



Eco-Friendly



Antimicrobial Surface



Hearing Aid Compatible



Distributed by:

SONIC
Everyday Sounds Better

1/10 Dorsey Street, Milton, QLD 4064
PO Box 1869, Toowong DC, QLD 4066
Free call: 1800 639 263 F 07 3368 1044
E enquiry@sonici.com.au
W www.sonici.com.au