

See It. Share It. Save It.



Welch Allyn Digital MacroView™ Otoscope

The Welch Allyn Digital MacroView Otoscope is a simple, affordable tool that lets you see, explain, and save digital images of the ear.

Offering you clear, live and still images of the ear canal and tympanic membrane, the Digital MacroView™ Otoscope is a valuable tool in support of the treatment and management of your patients' ears. Easy to learn and use, the Digital MacroView is a plug and play system requiring no additional purchase of software.

Here's how the Welch Allyn Digital MacroView Otoscope can benefit you:

- > Improves communication which can lead to better patient compliance
- > Allows you to store digital images in patient records for enhanced documentation of exam findings
- > Provides you the ability to share images with other providers as visual descriptive reference
- > Gives you an excellent training aid for educational purposes
- > Offers you the ability to view small details of tubes and hearing aids



WelchAllyn®

Advancing Frontline Care™

The Welch Allyn Digital MacroView™ Otoscope

Features you expect from an otoscope and more:



USING THE DIGITAL MACROVIEW OTOSCOPE IS EASY AS 1 – 2 – 3

1 Plug In.

Plug the otoscope into your computer using the USB cable.

2 Power On.

Put on ear specula and turn on light.

3 Play.

View live picture or capture image.

Ordering Information

23920 Digital MacroView Otoscope with 10 ft USB cable and set of four reusable ear specula

Accessories

23804	Insufflator Bulb	52434-U	4.25 mm Single use Ear Specula
24400-U	Additional set of four Reusable Specula	06500-U	Lamp
24302-U	2.5 mm Universal Poly Reusable Speculum	23900-10	Replacement USB Cable, 10 ft
24303-U	3.0 mm Universal Poly Reusable Speculum	71000-A	NiCad Rechargeable Handle (Plug-in)
23404-U	4.0 mm Universal Poly Reusable Speculum	71000-C	NiCad/Convertible Rechargeable Handle (Plug-in)
23405-U	5.0 mm Universal Poly Reusable Speculum	71670	NiCad Rechargeable Handle for Universal Charger
24320	SofSeal Small, Disposable Sealing Tip	71900	Lithium Ion Rechargeable Handle with AC Charging Module
24330	SofSeal Medium, Disposable Sealing Tip	71910	Lithium Ion Rechargeable Handle for Universal Charger
52432-U	2.75 mm Single use Ear Specula		

Distributed by:



1/10 Dorsey Street, Milton, QLD 4064
 PO Box 1869, Toowong DC, QLD 4066
Free call: 1800 639 263 F 07 3368 1044
E enquiry@sonici.com.au W www.sonici.com.au



Advancing Frontline Care™