



Workcover
State:
Claim No:

VIP ORDER
3 days in house

Custom Products Order Form

YOUR INFORMATION	Sonic Reference No.: <input type="text"/>			CLINICIAN'S INFORMATION		
	Phone/Email:			Today's Date:		
	Clinic Name:			Fitting Date:		
	Address:			Clinician's Name:		
				Clinician's Email:		
				CLIENT INFORMATION		
City:		State:	P/Code:	Full Name:		

CUSTOMS WITH FULL FLEXIBILITY	SELECT MODEL: Enchant¹⁰⁰ <input type="checkbox"/> Enchant⁸⁰ <input type="checkbox"/> Enchant⁶⁰ <input type="checkbox"/> Enchant⁴⁰ <input type="checkbox"/> Enchant²⁰ <input type="checkbox"/>								
	<input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> Both								
		Battery	NFMI (Binaural)	Wireless 2.4 GHz	Push button	Volume wheel	Telecoil	AutoPhone	Micro-phones
	<input type="checkbox"/> IIC	10	NA	NA	NA	NA	NA	NA	Omni
	<input type="checkbox"/> CIC	10	<input type="checkbox"/>	NA	<input type="checkbox"/>	NA	NA	NA	Omni
	<input type="checkbox"/> ITC	312	YES	<input type="checkbox"/> *	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> **	<input type="checkbox"/>	Directional
	<input type="checkbox"/> ITE HS	312 <input type="checkbox"/> 13 <input type="checkbox"/>	YES	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> **	<input type="checkbox"/>	Directional
<input type="checkbox"/> ITE FS	312 <input type="checkbox"/> 13 <input type="checkbox"/>	YES	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> **	<input type="checkbox"/>	Directional	
*ITC 312 2.4 GHz Depending on ear shape		<input type="checkbox"/> If ITC not possible do ITE with 2.4 GHz <input type="checkbox"/> If ITC not possible remove 2.4 GHz			** 2.4 GHz and Telecoil combination not available				

NFMI (Near-Field Magnetic Induction) enables a continuous exchange of data and audio between two hearing aids to provide advanced binaural processing

REQUIRED: AUDIOMETRIC INFORMATION										WIRELESS ACCESSORIES		
Hz	125	250	500	1K	2K	3K	4K	6K	8K	Hz	1K	
AC Right										BC Right		<input type="checkbox"/> SoundClip A
AC Left										BC Left		<input type="checkbox"/> Remote Control 3.0
												<input type="checkbox"/> TV adapter 3.0

FITTING LEVEL	75	85	90	100	SONIC SELECTS FITTING LEVEL FROM AUDIOGRAM (default) <input type="checkbox"/>		
IIC / CIC	<input type="checkbox"/>	<input type="checkbox"/>			WAXFILTER: Prowax MiniFit (default on all styles)		
ITC / ITE	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

SHELL COLOURS	BEIGE	MEDIUM BROWN	IIC COLOURS	RED/BLUE	CLEAR
FACEPLATE	<input type="checkbox"/>	<input type="checkbox"/>	FACEPLATE (Black)	<input type="checkbox"/>	<input type="checkbox"/>

VENTING Please select: OPTION A - OR - OPTION B																																												
A) LET SONIC SELECT		B) SPECIFY VENT																																										
R <input type="checkbox"/> L <input type="checkbox"/> Sonic Optimise Vent Size and Style (Audiogram Required) (Sonic will optimise venting to provide appropriate requirements for hearing loss and venting whilst minimizing instrument size)		<table border="1"> <thead> <tr> <th colspan="4">VENT EFFECT</th> <th colspan="2">VENT STYLE</th> </tr> <tr> <th>R</th><th>L</th><th>R</th><th>L</th> <th>R</th><th>L</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td><td><input type="checkbox"/></td> <td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> <tr> <td></td><td></td><td>0.8mm[^]</td><td>3.0mm</td> <td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> <tr> <td></td><td></td><td>1.4mm</td><td>4.0mm[^]</td> <td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> <tr> <td></td><td></td><td>1.8mm</td><td>No Vent</td> <td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> <tr> <td></td><td></td><td>2.4mm</td><td></td> <td><input type="checkbox"/></td><td><input type="checkbox"/></td> </tr> </tbody> </table> <p><small>^Size not available for collection vent (vent size can be reduced)</small></p>	VENT EFFECT				VENT STYLE		R	L	R	L	R	L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			0.8mm [^]	3.0mm	<input type="checkbox"/>	<input type="checkbox"/>			1.4mm	4.0mm [^]	<input type="checkbox"/>	<input type="checkbox"/>			1.8mm	No Vent	<input type="checkbox"/>	<input type="checkbox"/>			2.4mm		<input type="checkbox"/>	<input type="checkbox"/>
VENT EFFECT				VENT STYLE																																								
R	L	R	L	R	L																																							
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																							
		0.8mm [^]	3.0mm	<input type="checkbox"/>	<input type="checkbox"/>																																							
		1.4mm	4.0mm [^]	<input type="checkbox"/>	<input type="checkbox"/>																																							
		1.8mm	No Vent	<input type="checkbox"/>	<input type="checkbox"/>																																							
		2.4mm		<input type="checkbox"/>	<input type="checkbox"/>																																							

OPTIONS		
CANAL LENGTH (Factory standard unless otherwise marked)		EXTENDED WARRANTY (please see price list for charges)
<input type="checkbox"/> Medium <input type="checkbox"/> Long <input type="checkbox"/> As marked on impression		<input type="checkbox"/> Extra 1 year <input type="checkbox"/> Extra 2 years
OTHER OPTIONS		Note: All Sonic hearing aids come with a standard 3 Yr warranty (except HSP free-to-client models)
<input type="checkbox"/> Helix Lock	<input type="checkbox"/> Nail grip	
<input type="checkbox"/> Canal Lock	<input type="checkbox"/> Pull out string	
		VOLUME WHEEL
		<input type="checkbox"/> Standard^^ <input type="checkbox"/> Raised, standard^^ <input type="checkbox"/> Large^^ <input type="checkbox"/> Raised, large^^
		<small>^^Not available for size 10 instruments</small>

COMMENTS _____

SUBMIT FORM