



## WHAT IS HEARSCOPE?

A world-first, affordable video-otoscope for a smartphone. Take high-quality images and videos of an eardrum.

## ARTIFICIAL INTELLIGENCE

The user-friendly application means anyone, anywhere in the world, can record images on a mobile phone and upload for automated diagnostic support.

Our Artificial Intelligence system applies machine learning to provide an accurate diagnosis for the most common ear diseases.

\*Available soon



Reusable Specula



hearScope hardware with FREE app | Available on Android & USB computer connection

## KEY FEATURES

### VARIABLE LIGHT SOURCE

Light emitting diode (LED) bulbs.

### HIGH QUALITY 2MP CAMERA

Clear view for accurate diagnosis.

### 30 FRAMES/SEC VIDEO FRAME RATE

Real-time viewing and video capturing without delays.

### SELECTION OF SPECULA

Multiple specula included. Convenient range of specula to fit ears of all sizes (2x 3mm; 2x 4mm; 2x 5mm).

### SECURE CLOUD-BASED STORAGE

Sync hearScope images to the cloud to share with patients and colleagues.

\*Available soon

### VARIABLE MAGNIFICATION

Zoom like never before with your video-otoscope.

### DERMAL AND THROAT ADAPTOR

Attaches like a specula to take images of moles and skin lesions or fit tongue depressors to capture throat inspections.

+1 415 212 5500 | +27 (0)12 030 0268 | [info@hearXgroup.com](mailto:info@hearXgroup.com) | [www.hearXgroup.com](http://www.hearXgroup.com)

