

See The New MedRx Studio Software



One Great Software for all the MedRx Products


Simple to use and easy to customize for your clinic

MedRx[®]
www.medrx-usa.com

The New MedRx Studio Software

8 user-friendly modules


Incorporate the New *MedRx Studio Software* into your daily practice to deliver quality professional hearing healthcare.



DUAL CHANNEL AUDIOMETRY

Perform air conduction, bone conduction, masking and speech testing using the mouse and/or keyboard with ease.

This clean, modern design is easy to use and functions as your complete Audiological




REAL EAR MEASUREMENTS

IEC/ANSI standard described REM functions are provided for: REUR, REUG, REIG, REAG, REAR, REOR, REOG and RECD.

M STUDIO Software


LSM



LIVE SPEECH MAPPING

Live Speech Mapping (LSM) has quickly become a successful approach to hearing instrument fitting.


IMP



TYMPANOMETRY

This function provides clear graphic and numeric displays of tympanogram and acoustic reflex data.

HIT



HEARING INSTRUMENT TESTING

The IEC/ANSI standard tests are available in an automatic sequence or for select measurements.


VO



VIDEO OTOSCOPE

Capture, view and store digital images. Retrieve results easily for follow-up visits.


HLS
MHA



HLS/MHA

Easily simulate hearing loss and hearing aids for patients.

TI



TINNOMETER

Revolutionary Tinnitus Assessment.



The MedRx Studio Software is easy to use and equips you with excellent counseling tools.

MedRx[®]
www.medrx-usa.com

1200 Starkey Rd. #105
Largo, FL 33771 U.S.A.
888-392-1234 • 727-584-9600
FAX: 727-584-9602
Email: sales@medrx-usa.com
www.medrx-usa.com

Australian Distributor:
Sonic Innovations
1800 639 263
equipment@sonici.com

