

Otosure

Portable PC-based automatic audiometer



The **Amplivox Otosure** is a PC-based automatic audiometer that provides a pro-active and portable solution to satisfy audiometric screening requirements.

Technology

Incorporating digital technology, high quality standard audiometric connections and USB connectivity, the **Amplivox Otosure** provides a reliable and robust 'plug and play' solution that can be used with desktop PC's, laptops and, or PC tablets running the Windows® operating system. User definable test options provide a customisable solution and the **Otosure** can also interface to a number of leading OH database applications.

Pro-active

All tests can be reviewed in real time and on-screen comparisons can be made with previous results. The **Amplivox Otosure** also has a user definable recall function and extensive data analysis features that enable 'at risk' individuals and user defined groups to be identified and monitored.

Portability

Pocket sized and a weight of just 165g provide an extremely portable audiometer that can be used for both static and mobile requirements.

The **Amplivox Otosure** is a very cost effective solution for the provision, storage and management analysis of audiometric test results.

- **Pocket sized audiometer**
- **USB connectivity**
- **'Plug and play' functionality**
- **Test results analysis**
- **User definable recall facility**
- **Identification of 'at risk' personnel**
- **In-depth management analysis**



PC not supplied
Audiocups are optional



specialists in healthcare products

Otosure

Portable PC-based automatic audiometer



Technical specifications

Output/input data	
Frequency range (Hz):	125, 250, 500, 1K, 1K5, 2K, 3K, 4K, 6K, 8K
Frequency accuracy:	<1%
Distortion:	<2%
Output level range:	-10dBHL to 100dBHL (max)
Output level accuracy:	±3dB
Test modes:	Computer Manual
Output level step size:	Computer; 5dB Manual; 5dB
Threshold level:	Computer; 5dB Manual; 5dB

System requirements	
Operating system:	Stand alone PC with Windows XP or higher
Disc drives:	CD and hard disc drive
Memory:	Internal memory and available disc space as required by the PC operating system
Software:	The programme will be supplied on a CD
Interface:	USB

Physical data	
Power:	Via USB connection from PC
Dimensions (mm):	120 x 86 x 28 high
Weight:	165 grams

Standards	
The Amplivox Otosure conforms to the relevant clauses of the following standards.	
Safety:	IEC 60601-1
EMC:	IEC 60601-1-2
Audiometric:	IEC 60645-1
CE mark:	Complies to EU Medical Device Directive

Equipment

- | Standard equipment | |
|---|--|
| <ul style="list-style-type: none">● Audiometric headset● Patient response switch● Operation manual (CD)● USB cable● Custom-designed carrying case | |



Carrying case
For domiciliary use, the **Amplivox Otosure** is supplied with a carrying case, designed to accommodate the audiometer and all accessories.

- | Optional equipment | |
|---|--|
| <ul style="list-style-type: none">● Audiocups (noise reducing headset enclosures) | |

Australian Distributer:



SONIC
Everyday Sounds Better

Free Call:
1800 639 263

amplivox

Amplivox Ltd, 29-30 Station Approach
Kidlington, Oxford, United Kingdom OX5 1JD
Tel: +44 (0)1865 842411 • Fax: +44 (0)1865 841853
E-mail: sales@amplivox.ltd.uk

www.amplivox.ltd.uk